Laboratory Notebook - TPOB GPH4102

Laboratory Notebook - TPOB GPH4102

Raman Spectroscopy

Preparation

Theory

To do before starting

Necessary information

In Lab

Discussion with Simon Rainville

Manipulations

Software exploration

Readout Noise characterization

Thermal noise characterization

Photon noise (statistical)

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Second period

Substances Raman spectrum

Discussion with Daniel and Simon about Fluo/Raman

Confocal microscopy

Preparation

In Lab

Observations montage

Manipulations

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Résolution z

Fluorescence

Hi-Lo Microscopy

Preparation

Hi-Lo

Microscopie

In Lab

SEMAINE 1 MICROSCOPIE

Optical tweezers

Preparation

In Lab

Oximetry

Preparation

In Lab

Raman Spectroscopy

Preparation

Theory

Raman effect: re-emission of radiation from an excited state of a sample. [raman spectroscopy](#)

Rayleigh scattering: re-emission at same wavelength

Stokes: loss of energy in vibrational states excitation of the sample, re-emission has a bigger wavelength

Anti-Stokes: synchronous phonon (vibrations) gives a virtual higher state of energy, re-emission has smaller wavelength

Vibration of certain molecules or molecule bonds are relative to the bond. Thus, the absorption in the vibrational states will be different for different bonds. It is possible to characterize a material looking at its bonds.

Raman scattering is really weak. It is important to reduce/minimize instrumentation noise and other sources like thermal, external photons, etc.

To do before starting

1. Download LightSpec software (should be with the camera)

Necessary information

1. Mercury Spectrum for CCD calibration

[Data of Hg spectrum](#)

Important peaks for calibration will be between 633nm and lower (around 800nm since a $\Delta\lambda = 230$ is rare) since we excite our sample with a 633nm laser. According to [this site](#), Stoke's shift is much higher than anti-stokes, since anti-stokes requires the superposition of a photon while the electron is in a virtual state. Stoke's shift occurs more naturally, since the excitation will lead to an eigen energy value.

Interesting peaks with high relative intensity: 404.656nm(1800), 435.833nm(4000), 546.074nm(1100), 614.950nm(1000), 690.752(250), 708.190(250).

The spectrum on the CCD should start at around 630nm, so it will be calibrated with the 546nm and 614nm peaks, which will be highly distinguishable.

2. Pixis 100 CCD

[Explication théorique](#)

[Pixis100 Datasheet](#)

The CCD has 2 modes:

- High well sensitivity : Well depth of 300 ke- (typical), 250 ke- (min)
- High well capacity : Well Depth of 1 Me- (typical), 750 ke- (min)

The probability that a photon excites an electron is a function of the spectrum. A python code here makes the necessary corrections on data taken with the PIXIS100-B CCD.

[img](#)

Figure - Quantum efficiency of the 100PIXIX

Code to correct the measurements of the CCD-PIXIS100-B

The number of photons necessary to increment a pixel by 1 will depend on the mode selected. A higher capacity will require a higher number of electrons. The camera captor is coded in 16 bits: 65535 values. Thus, if we use the typical values:

- High well sensitivity : $\frac{300000}{65535} \approx 4.57771$ /electron per pixel incrementation.
- High well capacity : $\frac{1000000}{65535} \approx 15.25902$ /electron per pixel incrementation.

[QE graph linear approximation](#)

3. Static Photon Noise / Noise Relative to Measure

The Datasheet specifies static noise. No measure-dependent noise seems to be in effect here.

For the PIXIS 100-B, there are 2 acquisition speed, which affect the noise:

- 100kHz: 3e- rms (typical), 5e- rms (max)
- 2 MHz: 11e- rms (typical), 16e- rms (max)

Assuming the noise is a varying *sinusoidal* signal. The peak noise will rather be $\sqrt{2} \cdot U_{rms}$

The **noise relative to the measure** @-80°C is typically 0.001 electron/photon/second

It doesn't mention a dependency to acquisition speed/mode. So the measurement noise is given by this relation:

[Graph of relation](#)

4. Vibration Frequency

Here is a [list](#) of the usual bonds with emission wavelength associated

The Raman wavelength is obtained by: $\nu = \frac{1}{632.8nm} - \frac{1}{\lambda(x)}$. With $\lambda(x) = 700nm$, it gives 1517.07 cm^{-1} .

With this [table](#), it is possible to predict the emission wavelength of different molecules, as these tables reveal the emission of specific bonds. The values were taken from [this table](#), [this document](#)

5. Noise Evolution With Time Integration

Here is a [table](#) ready for the Time/Noise measurements.

6. Chlorophyll Spectrum

[Values taken on this site](#)

Absorption peaks : 430nm and 662nm

Emission peaks : 669nm

7. Olive Oil

Chemical compound of Olive Oil is a triglyceride with 3 fatty-acids.

$C_{55}H_{98}O_6$ is an example of chemical formula for a fat.

Olive Oil Fatty-acids are unsaturated. There should be presence of C=C bonds.

Also, in the glycerol, C-C, C-O, C=O, C-H, bonds will also be present.

This gives the expected range of emission of these bond

- 13.8% saturated (C-C)

- 70% mono-unsaturated (C=C)
- 13.2% poly-unsaturated (many C=C)

In Lab

Discussion with Simon Rainville

1. Difference between fluorescence and Raman scattering

Fluorescence is the excitation of an electron at a specific energy level, which then leads to the energy loss in complex non-radiative processes. Fluorescence excitation happens only at specific wavelengths.

Raman scattering has the same output: a bigger wavelength, thus a lower energy photon. The difference is in the excitation processes. The electron doesn't have to be excited at a specific level. It will rather be excited at a virtual level, and then release its virtual energy into vibrational modes of certain electron bonds. These energy variations vary according to the bond excited.

2. Three types of noise

- **Measure noise:** the action of measuring generates noise in the process. It is constant a each measure.
- **Thermal noise:** There is a probability that (because of the Boltzmann distribution of energy) it is possible that an electron will be excited. That probability is a function of time. So for example, There is 0.005% of chance that an electron is excited per second because the thermal energy at -80celcius. (not exact numbers)
- **Schottky noise:** Physical phenomenon's are usually described by Poisson's distribution. It is the case for photon emission of a sample. It is possible that the number of photon emitted are not totally the same for two measures. For Poisson's distribution, The standard deviation of the emission of a photon is the square root of number of photons detected. Thus, because we are interested in the ratio between the average value and the deviation, we shall do the ratio. So the Signal Noise Ratio (SNR) = $\frac{N}{\sqrt{N}} = \sqrt{N}$

Manipulations

Setup exploration

Important parts of the setup :

- Notch filter : Blocks the laser emission of 633nm that could come from rayleigh scattering or the laser himself. Important, light collected must arrive normal to the lens.
- Vertical Slit : manages the resolution. a thinner slit means less light, but a thinner bar of light, the goal is to make the slit of light the width of a pixel, so no resolution is lost and the maximum light enters the spectrometer.

Software exploration

WINSpec32 is the software to acquire spectrum with the PIXIS100 CCD. Here are the different parameters used during all the experimentation:

- The rate was set to 100K Hz
- The mode is High Capacity
- The gain was 1

- The ASCII output is on
- The accumulation of spectrum will probably vary, but is set to 1, will notice if changed
- The integration time will change

Picture of the setup dialog box:

To acquire, the Acquire button is used.

Readout Noise characterization

In WINSpec32, we asked the software to output a ASCII file. It saves .txt files with (spectrum number, x, y) convention

Integration time (ms)	Filename [*.txt]
100	bruit_mesure_100ms
75	bruit_meusre_75ms
50	bruit_meusre_50ms
25	bruit_meusre_25ms
10	bruit_meusre_10ms
5	bruit_meusre_5ms
1	bruit_meusre_1ms
0.1	bruit_meusre_100us
0.001	bruit_meusre_1us

The average noise is approximately 600/65535 per pixel.

INFO Binning doesn't divide by the number of pixel, it adds. So the limit of the value on the software is $65535 \times 100 = 650k$.

Daniel verified our configuration, He changed the Rate from 100kHz to 2MHz. The Readout noise measures were retaken with that configuration so the name of the files don't change.

Thermal noise characterization

Thermal noise should appear on long integration time.

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Observations : We noticed that long exposition increased the noise proportionally to the integration time. Also, there were weird spikes observable in the files. Even with all the lights off, the spikes are there. Maybe bleeding filter???

Photon noise (statistical)

We think that this noise should increased proportionally to \sqrt{N}

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Observations : General shape of the spectrum stays the same

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Second period

Substances Raman spectrum

Substance	Spectrum File	Integration time[s]	Spectrum accumulation	
Olive oil	olive_oil.txt			
Sunflower oil	sunflower_oil.txt	60	5	
Peanut Oil	peanut_oil.txt	10	10	
Canola Oil	canola_oil.txt	60	10	
Corn Oil	corn_oil.txt	60	5	
Ethanol	ethanol.txt	50	1	
Methanol	methanol.txt	180	1	
Isopropanol	isopropanol.txt	120	1	
Sucrose	sucrose.txt.	120	3	
Glycerol	glucose.txt	120	1	
A,B,C,D,E,F,G,H	A.txt, B.txt, C.txt, etc.	20	1	

For C accumulation was 3, for D 8, for G 5.

For all these substances, gain was set to 1. Rate was 2MHz. Readout was "Low Noise"

Discussion with Daniel and Simon about Fluo/Raman

Daniel explained to us a fundamental thing about measuring Raman spectrum in the presence of fluorescence. Since the signal of fluorescence is much more important than the Raman signal. However, it is possible to detect the Raman signal in the fluorescence signal. Since our Fluorescence signal is very large , for example the olive oil: the max signal after 100ms of integrating is 12000. The fluorescence signal is very strong. No peaks on the fluorescence signal can be seen.

If we take a pure sample, let's say ethanol. For 50s of integration, we got a minimum signal of about 2 000 as well. Because we want to see those small peaks

Thus, the signal is approximatively 3000 times weaker. for 100ms the Raman is 4. The statistical noise of the fluorescence is $\sqrt{12000} = 110$. We want a Raman signal approx. 5 times bigger than the noise.

$$4t[ms]/100 = 5 * \sqrt{12000 * t/100}$$

isolating t , it gives approx.. 30min of integration for a 5 times bigger Raman signal over the noise of the fluorescence. Retaking the spectrum of the olive oil effectively shows some spikes that weren't seen before.

It has been found that the high peaks at about 200pixels is the filter bleeding. The 4f system before the spectrometer was not aligned, which caused the light of the fluorescence to enter the system at a weird angle. The filter must have collimated light at its input

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